

PCR test kit information



Inspection equipment

CronoSTAR™ 96 Real-Time PCR System

A test-confirmed device capable of measuring new coronavirus reagents.



reagent

SARS-COV-2 Detection Kit Multi NCV-403

Real-time 1-step RT-PCR kit for detecting new coronavirus.



Saliva collection set

SG-915

A collection kit that can collect samples of the new coronavirus.

"Kinoshita Group New Corona PCR Testing Center" uses "Toyo Spinning Co., Ltd. New Corona Virus Detection Kit SARS- CoV

-2 Detection Kit" .. This product was announced by the Tuberculosis Infectious Diseases Division, Health Bureau, Ministry of Health, Labor and Welfare and the National Institute of Infectious Diseases as "2019nCoV genetic test method for which evaluation results using clinical specimens were obtained (September 30, 2020 version)". The result of obtaining the positive and negative concordance rates with the National Institute of Infectious Diseases is a product with both positive and negative concordance rates of 100%. This test is a genetic test using the real-time PCR method using this product.

- The tests conducted by this laboratory determine whether the new coronavirus is contained in saliva or nasal swab, and do not diagnose the morbidity of the new coronavirus infection.
- The purpose of the inspection service in stores is to reduce the burden on those who need inspection even though they do not need medical treatment, such as regular inspections as part of the prevention of infectious diseases. It is provided by. Therefore, for those who have symptoms such as fever and those who have close contact with us, we may refuse the examination and recommend that you visit a medical institution.
- This "sanitary laboratory" conducts PCR tests for asymptomatic persons, and if a positive (suspected) test is made, be sure to consult a medical institution for diagnosis. If you get a positive diagnosis at a medical institution, you will be notified to a health center, etc., so please seek medical attention immediately.